

MARS **HD**

Explore new world

Perfect for students and professionals, this versatile video magnifier offers 3 viewing modes making it easy to view objects close-up, materials on your desk, or presentations across the room!



 3 viewing modes

 HD camera

 Less than 0.9kg

ZOOMAX[®]
www.zoomax.com

Technical Specifications

Camera	12X optical zoom, autofocus HD with resolution 1280 x 720 pixels
Frame rate	60 full digital frame /second @ 720p
Magnification	2X ~ 80X (based on 14"screen)
Color modes	Full color + 16 high contrast color modes
3 viewing modes	Close, distance and self-viewing
Ergonomic adjustment	Rotate 330° vertically or 180° horizontally for multi-angles reading
USB/HDMI connection	Connect to Mac, PC and monitor
Functions	Zoom in/out, Color mode adjustment, Brightness adjustment, Reading lines/windows, Freeze frame, Split screen, Find, Image saving, Led light on/off, Setting memory
Battery	Around 4 hours continuous use, rechargeable Li-ion battery
Dimensions (folded)	22.6cm×8.6cm×12.5cm
Weight	Less than 0.9kg
Certificates	CE; FCC; RoHS



Contact us



Tel: +86-571-87006308

Email: sales@zoomax.com

Web: www.zoomax.com