

PORTABLE DEVICE

Snow Pad

Witness Every Moment



Snow Pad is the best option for visually impaired students to gain greater access in the classroom. It is also an ideal vision companion for every low-vision person to handle daily multi-tasks!

- OCR & Text-To-Speech
- Portable with Excellent Image Quality
- Compatible with Your Personal iPad



Technical Specifications

Magnification	1 – 32X
Display Modes	Full color and 10 high contrast color modes Auto
Build-in Distance Camera	Focus, Optical Zoom 7X, 1920 × 1080P, Full HD
Frame Rate	60FPS
Folded Dimensions	L 309*W 312*H 55mm
Stand Height	296mm
Weight	Approx. 2.4kgs (5.29lbs) without tablet
Battery	Approx. 5h
Charging Time	Approx. 4.5h
Power Adapter	Input 100 – 240V; Output 5V/3A
OCR/Text-to-speech	24 Languages
Compatible W/Tablets	All kinds of iPad; Xiaomi Pad 6 & Xiaomi Pad 6 Pro
Viewing Modes	4 Options: Split Screen; Distance-viewing; Near-viewing; Selfie

