

ZOOMAX[®]



LOW VISION AIDS



Zoom Your Vision, Max Your Life

About Us

Your Vision Product Expert

Established in 2011, ZOOMAX is a high-tech visual aid solution company. With long-term cooperation with famous universities such as the University of California (UCLA), Massachusetts Institute of Technology (MIT), etc., we are proud of our creativity in the hardware, software, and product ID search and design.

Why Choose Zoomax

- 14 years low vision industry experience
- Over 60 technology and product patents
- Certificate & Warranty: FDA ·RoHS ·CE ·FCC
- Distribution network in 80 countries worldwide
- Full line of low vision aids: handheld & portable, wearable & CCTV

Zoomax, User Experience Orientation

Zoomax pays high attention to the user's experience of the product: we incorporate user suggestions into the new product development. Zoomax people believe: only by loving and understanding our customers deeply, can we create customer-demanding products. Therefore, we always keep the dialogue with our clients and partners. Any valuable comments and suggestions are regarded as of great importance to us!

CONTENT

01 PORTABLE DEVICE 4-12

Snow Pad NEW	4
Snow Eye NEW	7
Snow 12	10

02 HANDHELD DEVICE 13-21

Luna S	13
Luna 6	16
Luna 8	19

03 DESKTOP DEVICE 22-27

Luna Eye NEW	22
Luna HD 24 Pro	25

04 WEARABLE DEVICE 28-30

Acesight VR	28
-------------------	----

PORTABLE DEVICE

Snow Pad

Witness Every Moment



Snow Pad is the best option for visually impaired students to gain greater access in the classroom. It is also an ideal vision companion for every low-vision person to handle daily multi-tasks!

- OCR & Text-To-Speech
- Portable with Excellent Image Quality
- Compatible with Your Personal iPad



01 Technical Specifications

Magnification	1 – 32X
Display Modes	Full color and 10 high contrast color modes
Build-in Distance Camera	Auto Focus, Optical Zoom 7X, 1920 × 1080P, Full HD
Frame Rate	60FPS
Folded Dimensions	L*W*H/ 309*312*55mm
Stand Height	278mm
Weight	Approx. 2.4kgs (5.29lbs) without tablet
Battery	Approx. 5.5h
Charging Time	Approx. 4.5h
Power Adapter	Input 100 – 240V; Output 5V/2.5A
OCR/Text-to-speech	20 Languages
Compatible W/Tablets	All kinds of iPad; Xiaomi Pad 6 & Xiaomi Pad 6 Pro
Viewing Modes	4 Options: Split Screen; Distance-viewing; Near-viewing; Selfie



02 Users' Reviews

I was happy to have the opportunity to test the Snow Pad some time ago. Firstly, I love the compact design and the high-quality camera. My patients tried it out and expressed their satisfaction for the easy-to-use features. I believe many low-vision people who have iPad will deeply love Snow Pad!

Doctor – April. 2024, USA

Snow Pad was the most impressive product I found at CSUN 2024 conference. The distance-viewing is so important for the students with visual impairments. The compatibility with iPad is an amazing idea that helps take great advantage of the students' personal tablet. I had recommended Snow Pad to my students and other teachers.

Teacher – March, 2024 USA



PORTABLE DEVICE

Snow Eye Eye on your love



- Connect with Snow 12 or any screens with HDMI
- Plug and start your distance-view instantly



Snow Eye



01 Technical Specifications

Magnification	1X – 32X
Display Modes	Full color and 12 high contrast color modes
Camera	Optical Zoom, Auto Focus, 1920 × 1080p, 60FPS
Weight	350g (0.77lbs)
Dimension	92*92*117mm (3.6*3.6*4.6in)
Battery	Approx 4.5hrs continuous use
Charging Time	Approx 3.5hrs
Connection	Connect with monitor, laptop and Zoomax Snow 12



02 Users' Reviews

We tested the Snow Eye distance camera in our center. Many users found it a great help in daily life. Some of them connected Snow Eye with their computer, while some of them use laptop or Zoomax Snow 12. It doesn't need any firmware installation and definitely easy for the setup. The optical zoom brings amazing picture quality! We all love it!

Low Vision Center – April, 2024, USA

Snow Eye is my most favorite low vision aids at Sightcity 2024! I love its compact design and extremely superior image quality. My vision is very bad. It is unbelievable that I could see things ten meters away clearly with the help of Snow Eye! Really hope it will benefit more people with vision loss.

Thomas – May, 2024, Germany



PORTABLE DEVICE

Snow 12 Simple and Powerful



As a portable electronic video magnifier, Snow 12 features a simple accessible interface, which makes it easier for users to read, write or conduct full-page content tasks.

- OCR & Text-To-Speech
- Detachable & Foldable for Outdoor Use
- Button or Touchscreen with Voice Prompt



01 Technical Specifications

Magnification	2.7X -19X; 2.5X -19X w/stand
Color Modes	Full color & 10 high contrast color modes
OCR and Text-to-Speech (optional)	More than 20 languages
File Formats	JPG, PDF, RTF, and TXT
Screen	12-inch, 1920*1080
Camera	13M pixels
Weight	Snow 12: 960g (2.1lbs); 1230g (2.7lbs) w/stand
Dimensions	29*22*2.4 cm (12*9*1 in.)
Continuous Use Time	2.5 to 3 hours
Charging Time	2 hours
Storage Capacity	13 GB



02 Users' Reviews

I bought Snow 12 for my grandma. She said it performs very well. Whenever she had any questions, she talked with the Zoomax people and they always help her in time.

Jonathan - March 4, 2022, USA

My vision is very limited and Snow 12 is a life-saver. It allows me to read the daily newspaper without holding a magnifying glass. I leave it on my dining table for convenience. I also take it when I go out, so I can use it anytime I like.

Lloyd - July 10, 2022, New Zealand



HANDHELD DEVICE

Luna S

Better and Happier



The 4.3-inch Luna S is a useful tool for those with low vision. Luna S offers sharp photos and outstanding features – it can be used to read menus, books, magazines, tags, and even listen to the radio.

- Crystal-clear Images
- FM Radio in the Pocket
- Portable and Lightweight



01

Technical Specifications

Magnification	2X - 19X
Display Modes	Full color and 10 high contrast color modes
Focus	Auto focus
Screen	4.3 inch
Flashlight Mode	Yes
FM Radio	Yes
Voice Control	Available in specific countries and regions
Power Saving	Unit will turn off automatically after no activity for 3 minutes
Battery	Replaceable lithium-ion battery
Battery Life	3 hours
Charging Time	4 hours
Weight	200g / 7oz
Dimensions	16cm x 8cm x 2.4cm / 6.3" x 3.15" x 0.94"



02 Users' Reviews

I'm surprised by the FM radio function of Luna S. It's not only the best visual aid I've ever used but also a product of entertainment for me. I always use Luna S to listen to my favorite radio station. And the matching free ear bugs are of high quality.

Tracy J - March 23, 2022, USA

Luna S is definitely the lightest, the most convenient hand-held electronic magnifier I've ever used. And Zoomax people are all responsible and friendly!

Mia - March 29, 2022, Singapore

My brother needs a visual aid at school. The customer service recommended Luna S, and we were so satisfied with the image and the design. Thanks for your great recommendation!

Wayne - August 8, 2022, UK



HANDHELD DEVICE

Luna 6

We Design It as You Need



Pocket-sized Luna 6 has a unique wireless charging system and built-in handle, offers sharp images even under fast movement. Users can easily carry it around and use it everywhere.

- Pocket-size and Lightweight
- Customized Wireless Charger
- Comfortable Built-in Handle



01 Technical Specifications

Magnification	2.5X - 19X
Screen	5.5-inch, Resolution 1280*720
Camera	13MP
Focus	Auto focus and refocus by short pressing the Mode button and Freeze button simultaneously
Color Modes	Full color and 10 high contrast color modes
Weight	Luna 6: 210g/0.46lbs; wireless charger: 92g/0.21lbs
Dimensions	160 x 80 x 11mm/6.3 x 3.1 x 0.4in
Battery Life	4 hours continuous use
Wired Charging Time	Less than 3 hours
Wireless Charging Time	About 4 hours




02 Users' Reviews

I bought Luna 6 from the clinic, and the service staff helped me a lot with the operation. I don't know how to express my gratitude to Zoomax. Thanks for making such a good product.

Tamara - April 5, 2022, Britain

I work in the library. I always need to check the book label. Luna 6 is the best device for me! It is light and slim so that I can bring it and use it anytime.

Kate - June 23, 2022, USA



HANDHELD DEVICE

Luna 8

Simple, Useful, and Stylish



Handheld electronic video magnifier Luna 8 is ideal for using anywhere owing to its small size, lightweight design, and 8 hours of battery life. It is friendly for both young and elderly low vision users.

- Sharp Image Instantly
- 8 Hours Long Battery Life
- Easy Operation & Use Everywhere



Luna 8



01 Technical Specifications

Magnification	2.5X - 19X
Color Modes	Full color and 10 high contrast color modes
Screen	8-inch, Resolution 1280*800
Camera	Auto focus; 13MP
Weight	500g (1.1lbs)
Dimensions	218 x 143 x 11mm / 8.6 x 5.6 x 0.4in
Battery Life	8 hours
Charging Time	Less than 3 hours



02 Users' Reviews

Luna 8 is absolutely helpful for my Mom. The image is super clear. It's definitely a nice-looking item. My Mom likes it a lot!

Christen - September 27, 2022, USA

The battery life is fantastic! I used Luna 8 once a day for half an hour. It had been working for almost 15 days with fully charge.

Margret - November 21, 2022, Britain

Luna 8 looks like a cool tech gadget and it's easy to use. It has a distance viewing function, so I can see the menu board in a restaurant or the slides in a conference room easily!

Donal - August 22, 2022, USA



CCTV DESKTOP DEVICE

Luna Eye

For All Your Vision Needs



Luna Eye is lightweight with excellent flexibility. You will easily lift it up and move it freely in your house or office. Seeing things far away, enjoying your favorite book, applying make-up or personal grooming - everything is possible!

- Distance and Near View; Self-view
- Easy-To-Use and Be Transported Easily
- Fast Auto-Focus and Superb Image



01 Technical Specifications

Camera	11x optical zoom lens, auto focus, 1920*1080, 60FPS
Magnification Range	1x - 72x
Screen	24-inch monitor with resolution of 1920*1080
Color Mode	True color + 18 kinds of high contrast color mode
Dimension	LWH 540*302*626mm, writing height Approx. 520mm
Basic Functions	Zoom in/out, color mode adjustment, panning, freeze
Other Function	Find function, reading line and mask, focus lock, contrast adjustment, brightness adjustment
Weight	approx. 6.8kgs/14.9lbs
Viewing Modes	Distance-view, near-view, selfie




02 Users' Reviews

My grandma has been using Zoomax Snow for ten years. Luna Eye is her new purchase. She reads and writes everyday. She is happy with the larger-size screen. She said the distance-viewing feature is also cool when looking for some tiny things around her.

Arun – May, 2024, India

We have been looking for an all-in-one system that can be used by our users with various daily needs. Zoomax Luna Eye is absolutely one of the best options. Not bulky, but compact and lightweight, and sharper image, Luna Eye is the key low vision product we will recommend to our users this year.

Low Vision Center – May, 2024, USA



CCTV DESKTOP DEVICE

Luna HD 24 Pro

Enhance Your Vision Experience



Luna HD 24 Pro is a stable, high-quality, and powerful full HD desktop video magnifier. Luna HD 24 Pro provides remarkable clarity high contrast and brilliant image, it's the best solution for users to enjoy reading, writing, and crafting.

- 3 Tactile Buttons Only, Plug & Play
- Excellent Image and 24-Inch Full HD
- Adjustable Bracket for Ultimate Comfort



Luna HD 24 Pro



01 Technical Specifications

Magnification	2.5x - 70x
Color Mode	True color + 20 kinds of color mode
Writing Height	Approx. 17.5cm
Camera	Optical zoom lens, auto-focus, 1920*1080P
Screen	24 - inch monitor with resolution of 1920*1080
Weight	Approx. 15Kg
Basic Function	Zoom in/out, color mode adjustment, find function, contrast adjustment, brightness adjustment, etc.
Advanced Function	20+ colors modes adjustment, find feature, reading line and mask , focus lock, freeze frame, memory setting, etc.



02 Users' Reviews

I can't believe that there is such a great product for low vision. With the help of Luna HD 24 Pro, I am able to continue the work I love and restart my hobby as well. I just can't stop expressing my heartfelt thanks to you!

Jasper - October 30, 2022, Russia

After trying Luna HD 24 Pro, decided to buy it. I like the large-screen electronic magnifiers, and Luna HD 24 Pro is the most user-friendly one among them! Even at large magnification, the images and text are still clear.

Sara - November 19, 2022, United Arab Emirates



WEARABLE DEVICE

Acesight 



As all-in-one ease-to-use wearable smart glasses, Acesight VR is a fantastic option for staying connected - advanced visual technology ensures better reading and watching experiences.



Acesight VR

Suitable for the following conditions:

- ✓ Macular Degeneration
- ✓ Glaucoma
- ✓ Diabetic Retinopathy
- ✓ Stargardt's Disease
- ✓ Nystagmus
- ✓ Esotropia
- ✓ Retinitis Pigmentosa
- ✓ And more

01 Technical Specifications

Magnification	1X - 16X
Camera	48 MP
Screen Resolution	Full HD 2340 x1080
FOV	65°
Refresh Rate	60 fps
Weight of Headset	460g/1 pound
Weight of Controller	65g/0.14 pound
Dimensions of Headset	18 x 12 x 10cm/7x4.7x3.9in
Headset Charging Time	Less than 1.5 hours
Headset Battery Life	2.5-3 hours



02 Users' Reviews

Acesight VR is so crucial for my husband. He is a handicraftsman and he needs to stare at small details sometimes. Acesight VR helps him a lot so he can continue his work. We cannot thank you enough!

Dennis - March 2, 2022, USA

It's a fantastic gift for my Dad. He likes Acesight VR! It's perfect for him to wear it while watching TV shows. That was what he dreamed of. I am so glad that I selected the right thing for him.

Dorothy - November 13, 2022, Canada





CONTACT US

Zoomax Technology Co., Limited

 info@zoomax.com

 www.zoomax.com

Zoomax (USA) Inc.

 info@zoomaxusa.com

 www.zoomaxusa.com