

# WEARABLE DEVICE

# Acesight



As all-in-one ease-to-use wearable smart glasses, Acesight VR is a fantastic option for staying connected – advanced visual technology ensures better reading and watching experiences.



Acesight VR

Suitable for the following conditions:

- |                        |                        |
|------------------------|------------------------|
| ✓ Macular Degeneration | ✓ Nystagmus            |
| ✓ Glaucoma             | ✓ Esotropia            |
| ✓ Diabetic Retinopathy | ✓ Retinitis Pigmentosa |
| ✓ Stargardt's Disease  | ✓ And more             |

# 01 Technical Specifications

|                       |                           |
|-----------------------|---------------------------|
| Magnification         | 1X - 16X                  |
| Camera                | 48MP                      |
| Screen Resolution     | Full HD 2340x1080         |
| FOV                   | 65°                       |
| Refresh Rate          | 60FPS                     |
| Weight of Headset     | 455g/1 pound              |
| Weight of Controller  | 65g/0.14 pound            |
| Dimensions of Headset | 180x120x100mm/7x4.7x3.9in |
| Headset Charging Time | Less than 1.5 hours       |
| Headset Battery Life  | 3.5-4 hours               |

# Users' Reviews

Acesight VR is so crucial for my husband. He is a handicraftsman and he needs to stare at small details sometimes. Acesight VR helps him a lot so he can continue his work. We cannot thank you enough!

Dennis - March 2, 2022, USA

It's a fantastic gift for my Dad. He likes Acesight VR! It's perfect for him to wear it while watching TV shows. That was what he dreamed of. I am so glad that I selected the right thing for him.

Dorothy - November 13, 2022, Canada

