

HANDHELD DEVICE

Luna 8

Simple, Useful, and Stylish



Handheld electronic video magnifier Luna 8 is ideal for using anywhere owing to its small size, lightweight design, and 8 hours of battery life. It is friendly for both young and elderly low vision users.

- Sharp Image Instantly
- 8 Hours Long Battery Life
- Easy Operation & Use Everywhere



Luna 8

01 Technical Specifications

Magnification	2.5X - 19X
Color Modes	Full color and 10 high contrast color modes
Screen	8-inch, Resolution 1280*800
Camera	Auto focus; 13MP
Weight	500g (1.1lbs)
Dimensions	218x143x11mm / 8.6x5.6x0.4in
Battery Life	8 hours
Charging Time	Less than 3 hours


02 Users' Reviews

Luna 8 is absolutely helpful for my Mom. The image is super clear. It's definitely a nice-looking item. My Mom likes it a lot!

Christen - September 27, 2022, USA

The battery life is fantastic! I used Luna 8 once a day for half an hour. It had been working for almost 15 days with fully charge.

Margret - November 21, 2022, Britain



Luna 8 looks like a cool tech gadget and it's easy to use. It has a distance viewing function, so I can see the menu board in a restaurant or the slides in a conference room easily!

Donal - August 22, 2022, USA