


# CCTV DESKTOP DEVICE

## Luna Eye

For All Your Vision Needs



Luna Eye is a 24" screen desktop electronic magnifier combining flexibility and space-saving advantage. Its compact, plug-and-play design ensures effortless use, while the intuitive interface enhances accessibility.

- Distance, near, and self-viewing
- Fast Auto-Focus and Superb Image



# 01 Technical Specifications

Camera	11X optical zoom lens, auto focus, 1920*1080, 60FPS
Magnification Range	1X - 72X
Screen	24-inch monitor with resolution of 1920*1080
Color Mode	True color + 18 kinds of high contrast color mode
Dimension	LWH 540*302*626mm, writing height Approx. 520mm
Basic Functions	Zoom in/out, color mode adjustment, panning, freeze
Other Function	Find function, reading line and mask, focus lock, contrast adjustment, brightness adjustment
Weight	approx. 6.8kgs/14.9lbs
Viewing Modes	Distance-view, near-view, selfie

## 02 Users' Reviews

Luna Eye is my grandma's new purchase. She reads and writes everyday, She is happy with the larger-size screen. She said the distance-viewing feature is also cool when looking for some tiny things around her.

Rani - January 15, 2025, India

We have been looking for an all-in-one system that can be used by our users with various daily needs. Zoomax Luna Eye is absolutely one of the best options. Not bulky, but compact and lightweight, and sharper image, Luna Eye is the key low vision product we will recommend to our users this year.

Low Vision Center - May 6, 2024, USA