

HANDHELD DEVICE

Luna S Better and Happier



The 4.3-inch Luna S is a useful tool for those with low vision. Luna S offers sharp photos and outstanding features – it can be used to read menus, books, magazines, tags, and even listen to the radio.

- Crystal-clear Images
 - FM Radio in the Pocket
- Portable and Lightweight



Luna S

01 Technical Specifications

Magnification	2X - 19X
Display Modes	Full color and 10 high contrast color modes
Focus	Auto focus
Screen	4.3 inch
Flashlight Mode	Yes
FM Radio	Yes
Voice Control	Available in specific countries and regions
Power Saving	Unit will turn off automatically after no activity for 3 minutes
Battery	Replaceable lithium-ion battery
Battery Life	3 hours
Charging Time	4 hours
Weight	200g / 7oz
Dimensions	160x80x24mm / 6.3x3.5x0.94in

02 Users' Reviews

I'm surprised by the FM radio function of Luna S. It's not only the best visual aid I've ever used but also a product of entertainment for me. I always use Luna S to listen to my favorite radio station. And the matching free ear bugs are of high quality.

Tracy J - March 23, 2022, USA

Luna S is definitely the lightest, the most convenient hand-held electronic magnifier I've ever used. And Zoomax people are all responsible and friendly!

Mia - March 29, 2022, Singapore

My brother needs a visual aid at school. The customer service recommended Luna S, and we were so satisfied with the image and the design. Thanks for your great recommendation!

Wayne - August 8, 2022, UK