

PORTABLE DEVICE

ZOOMAX[®]

Snow 12

Simple and Powerful



As a portable electronic video magnifier, Snow 12 features a simple accessible interface, which makes it easier for users to read, write or conduct full-page content tasks.

- OCR & Text-To-Speech
- Detachable & Foldable for Outdoor Use
- Button or Touchscreen with Voice Prompt



Snow 12

01 Technical Specifications

Magnification	2.7X - 19X/2.5X - 19X on the Foldable Stand
Display Modes	Full color and 10 high contrast color modes
Screen	12" touch screen
Screen Resolution	1920x1080
Panning	Touch screen panning & analog stick panning
Camera	13MP
Focus	Auto focus/tap to focus
Reading Line and Masks	Adjustable horizontal/vertical reading line and masks
Input Port	Type-C/HDMI
Battery	2.5-3 hours
Charging Time	2 hours
Standby Mode	Short press the Power button or stop using after 3 mins
PC Connection	Via Type-C port
TV/Monitor Connection	Via HDMI-in port

02 Users' Reviews

I bought Snow 12 for my grandma. She said it performs very well. Whenever she had any questions, she talked with the Zoomax people and they always help her in time.

Jonathan - March 4, 2022, USA

My vision is very limited and Snow 12 is a life-saver. It allows me to read the daily newspaper without holding a magnifying glass. I leave it on my dining table for convenience. I also take it when I go out, so I can use it anytime I like.

Lloyd - July 10, 2022, New Zealand

