

# PORTABLE DEVICE

## Snow Eye Eye on your love



Zoomax Snow Eye, a portable electronic magnifier designed to enhance distance viewing. Snow Eye connects easily to Snow 12, personal laptops or computers, enabling users to enjoy activities like watching TV, participating in conferences, or appreciating outdoor scenery.

- Connect with Snow 12 or any screens with HDMI
- Plug and start your distance-view instantly



Snow Eye

# 01

# Technical Specifications

|               |                                                 |
|---------------|-------------------------------------------------|
| Magnification | 1X – 32X                                        |
| Display Modes | Full color and 12 high contrast color modes     |
| Camera        | Optical Zoom, Auto Focus, 1920×1080, 60FPS      |
| Weight        | 350g ( 0.77lbs )                                |
| Dimension     | 92*92*117mm (3.6*3.6*4.6in)                     |
| Battery       | Approx. 4.5 hours continuous use                |
| Charging Time | Approx. 3.5 hours                               |
| Connection    | Connect with monitor, laptop and Zoomax Snow 12 |

# Users' Reviews

We tested the Snow Eye distance camera in our center. Many users found it a great help in daily life. Some of them connected Snow Eye with their computer, while some of them use laptop or Zoomax Snow 12. It doesn't need any firmware installation and definitely easy for the setup. The optical zoom brings amazing picture quality! We all love it!

Low Vision Center – April, 2024, USA

Snow Eye is my most favorite low vision aids at Sighticiy 2024! I love its compact design and extremely superior image quality. My vision is very bad. It is unbelievable that I could see things ten meters away clearly with the help of Snow Eye! Really hope it will benefit more people with vision loss.

Thomas – May, 2024, Germany

