

# Snow Pad 2025

Our latest foldable and portable video magnifier. The split screens with near and distant viewing mode empowers low vision individuals to handle multi-tasks. Outstanding compatibility allows for seamless integration with your personal tablets.



**Ideal For Home, School, Work & Travel**  
**Compatible With Your Personal Tablet**

**ZOOMAX<sup>®</sup>**



Introducing the Zoomax Snow Pad, a next-generation portable digital magnifier that seamlessly integrates with your personal tablet. Designed for individuals with low vision, the Snow Pad transforms your tablet into a versatile visual aid, offering exceptional image quality for both distance and close-up viewing. Its features include split-screen functionality, OCR/text-to-speech capabilities, and a customized app tailored for ease of use. Whether you're at school, work, or enjoying leisure activities, the Snow Pad enhances your visual experience, promoting independence and accessibility in daily life.



Customized App



OCR/Text to Speech



File Edit



Share File



Split Screen



Magnification

- ①Download Snow Pad App to Your Tablet
- ②Connect Snow Pad with Your Tablet using The Supplied USB Cable
- ③Power On The Snow Pad and Tap The Snow Pad App to Enjoy It

### Where to download the App?



Snow Pad iPadOS



Snow Pad Android



Turn 270 degrees clockwise.  
Turn 60 degrees counterclockwise.



Up 15 degrees.  
Down 40 degrees.

Magnification	1X – 32X
Display Modes	Full color and 10 high contrast color modes
Build-in Distance Camera	Auto Focus, Optical Zoom 7X, 1920×1080, Full HD
Frame Rate	60FPS
Folded Dimensions	L*W*H/ 309*312*55mm
Stand Height	278mm
Weight	Approx. 2.4kgs (5.29lbs) without tablet
Battery	Approx. 5h
Charging Time	Approx. 4.5h
Power Adapter	Input 100–240V; Output 5V/2.5A
OCR/Text-to-speech	20+ Languages
Compatible W/Tablets	iPad: iOS 16 or higher Android tablet: Xiaomi Pad 6, Xiaomi Pad 6 Pro; Galaxy Tab S9 FE, Galaxy Tab S10 FE
Viewing Modes	4 Options: Split Screen; Distance-viewing; Near-viewing; Selfie



**[www.zoomax.com](http://www.zoomax.com)**  
**[www.zoomaxusa.com](http://www.zoomaxusa.com)**  
**[info@zoomax.com](mailto:info@zoomax.com)**

