

# HANDHELD DEVICE

## Luna 6

### We Design It as You Need



Pocket-sized Luna 6 has a unique wireless charging system and built-in handle, offers sharp images even under fast movement. Users can easily carry it around and use it everywhere.

- Pocket-size and Lightweight
- Customized Wireless Charger
- Comfortable Built-in Handle



Luna 6

# 01 Technical Specifications

|                        |                                                                                           |
|------------------------|-------------------------------------------------------------------------------------------|
| Magnification          | 2.5X - 19X                                                                                |
| Screen                 | 5.5-inch, Resolution 1280*720                                                             |
| Camera                 | 13MP                                                                                      |
| Focus                  | Auto focus and refocus by short pressing the Mode button and Freeze button simultaneously |
| Color Modes            | Full color and 10 high contrast color modes                                               |
| Weight                 | Luna 6: 228g/0.5lbs; wireless charger: 92g/0.21lbs                                        |
| Dimensions             | 160x80x11mm/6.3x3.1x0.4in                                                                 |
| Battery Life           | 4 hours continuous use                                                                    |
| Wired Charging Time    | Less than 3 hours                                                                         |
| Wireless Charging Time | About 4 hours                                                                             |

## 02 Users' Reviews

I bought Luna 6 from the clinic, and the service staff helped me a lot with the operation. I don't know how to express my gratitude to Zoomax. Thanks for making such a good product.

Tamara - April 5, 2022, Britain

I work in the library. I always need to check the book label. Luna 6 is the best device for me! It is light and slim so that I can bring it and use it anytime.

Kate - June 23, 2022, USA

