

CCTV DESKTOP DEVICE

Luna HD 24 Pro

Enhance Your Vision Experience



Luna HD 24 Pro is a stable, high-quality, and powerful full HD desktop video magnifier. Luna HD 24 Pro provides remarkable clarity high contrast and brilliant image, it's the best solution for users to enjoy reading, writing, and crafting.

- 3 Tactile Buttons Only, Plug & Play
- Excellent Image and 24-Inch Full HD
- Adjustable Bracket for Ultimate Comfort



Luna HD 24 Pro

01 Technical Specifications

Magnification	2.5X - 70X
Color Mode	True color + 20 kinds of color mode
Writing Height	Approx. 17.5cm
Camera	Optical zoom lens, auto-focus, 1920*1080
Screen	24-inch monitor with resolution of 1920*1080
Weight	Approx. 15Kg
Basic Function	Zoom in/out, color mode adjustment, find function, contrast adjustment, brightness adjustment, etc.
Advanced Function	20+ colors modes adjustment, find feature, reading line and mask , focus lock, freeze frame, memory setting, etc.

02 Users' Reviews

I can't believe that there is such a great product for low vision. With the help of Luna HD 24 Pro, I am able to continue the work I love and restart my hobby as well. I just can't stop expressing my heartfelt thanks to you!

Jasper - October 30, 2022, Russia

After trying Luna HD 24 Pro, decided to buy it. I like the large-screen electronic magnifiers, and Luna HD 24 Pro is the most user-friendly one among them! Even at large magnification, the images and text are still clear.

Sara - November 19, 2022,
United Arab Emirates

